

Common Customisations & Accessories

Index

TA iQ Keep it black or spice it up	5
Rims & Tyres	5
Power Elevate, Tilt & Back Recline	6
Centre Legrests	7
Footplates	8
Individual Legrests	9
Amputation Legrests	10
Headrests	11
Leg supports/covers	11
Armrests	12
Joystick Controllers	13
Joystick Knobs	13
Alternative Controls	14
Attendant Controls	15
Flight Kit	16
LED Lights & USB Charger	16
Storage Bags & Hooks	17
Push Handles	17
Side Mirror	17
Cup & Device holders	18
Trays	19
Seat Pan Systems	20
Electric Seat Glide Function	20
Back Canes	21
Seat Belts	21
Chair Width	21
TA iQ Junior & TA Junior, Seating set-up accessories	22
Vehicle Docking System	23

We aim for a perfect fit every time...

We can upgrade your chair easily with many accessories and customisations.

This to ensure your daily tasks can be performed safely in the most comfortable way.

This brochure showcases the many adjustments and enhancements that can be implemented with TA's accessories. Morton & Perry technicians can offer further customisations on request.











Standard features at no upcharge included in every chair:



Seatplate to floor only 38cm

- > With electric HiLo elevate function. Standard on all chairs
- > Easier access under table and desks
- > Lower centre of gravity improves safety & stability at speed



Expandable Electronics

- > Standard access to 4-6 power functions capability see power functions on pages 6, 7, 8, 14, 15, 18
- > Access to a alternative controls capability page 13



Travel further

> No range anxiety when out and about TA iQ Models: 40km, TA Models 22km



Travel Safer

Occupant can remain in chair in vehicle

> Crash Test Safety Compliant - to ISO7176-19. NZ LVVTA 45-60 2.6



Full support from our team

- > Individualised sequencing programing
- > Simpler memory function setup
- > After sale support

TAM Keep it black or spice it up...









If you would like to add your own flavour;
Custom colours or designs are available upon request.

For the TA iQ models only. Additional charges my apply.

Rims & Tyres



TA iQ comes standard with black rims and tyres.

You can swap them out for alloy rims and/or swap the tyres out for gray.

(Custom colours can only be added to the black rims)







TA Indoor comes standard with alloy rims and grey tyres.



Power Elevate, Tilt & Back Recline

Included standard

Telescopic HiLo elevate function

The HiLo function is included standard with each TA & TAiQ chairs.

- It allows for a 30cm lift (20cm for the TAiQ standup)
- > Allows you to sit lower, so the chair fits under tables and desks. Also gives you more head clearance if you remain in the chair in a vehicle.




Standard with each TA & TAiQ* chair. Choose between manual or electric. The electric recline offers extra user comfort. Or Backrest recline Backrest recline code: XTA80363 Included standard Backrest power necline code: XTA97233



* The TA Standup comes standard with an electric set tilt of 45°, electric seat lift of 20cm, a recline of 180° lie flat and a centre electric legrest.

Centre Legrests



Both the Electric Centre Standard and Short Legrest have multiple adjustments available for length and angle.



Electric centre standard

Min Adjustment


Seat to Footplate up 13"/33cm Seat to Footplate down 18.5"/47cm

Max Adjustment

Seat to Footplate up 16"/41cm Seat to Footplate down 22"/56cm

Legrest power centre - std (with calf pads)

code: XTA80357



Electric centre short

Min Adjustment


Seat to Footplate up 10"/25cm Seat to Footplate down 13.5"/34cm

Max Adjustment

Seat to Footplate up 12.5"/32cm Seat to Footplate down 15.5"/39cm

Legrest power centre short with calf pads

code: XTA80571





Footplates

Standard



Footplate - aluminum with cut out edge (standard for MWD & RWD) code: XTA80540



Footplate - Standard complete (standard for Indoor & FWD) code: XTA80361-1

Special Footplates



Custom footplate - wider with extended lip

code: XTA80540-1



Footplate - Divided left & right pair.

code: XTA80818



Neoprene Padded Footplate code: XTACOVERS05

Customisations of footplates are available contact us to discuss your requirements.

Individual Legrests Amputation Legrests



Legrest - manual - RH s/a Stepless

code: XTA80700-1

Legrest - manual - LH s/a Stepless

code: XTA80701-1

Measurement from seat pan to foot plate


Shortest 29cm. Longest 50cm.

(It is possible to make longer tubes if needed. Custom sized legrest length available on request.)

Footplate and legrests adjustable to all angles.







 $\label{lem:lemma$


code: XTA80628

Legrest - electric elevating - LH with calf pad

code: XTA80629

Measurement from seat pan to foot plate

Shortest 22cm. Longest 50cm.





Headrests Leg support / covers



TA Standard Headrest Headrest system, cushion + bracket - Med 7" (18cm) code: XTA80337-1



TA Standard Headrest Headrest system, cushion, bracket STD 9" (23cm) code: XTA80337



Thigh support cushion (ea)
Med (10 x 13cm) code: XTA93263
Lrg (13 x 15cm) code: XTA93264
Custom sizes available on request



Nxt Head Support Contoured **code: XDNIB0785** 26x14cm



Nxt Head Support Large Rectangle code: XDNIB0784 28x14cm



Nxt Head Support XL code: XDNIB0818 38x14cm



Drop in seat guides code: XTA94000
Custom sizes available on request



TA centre electric neoprene cover (each) code: XCOVERS02



Swing away Legrest Neoprene Cover code: XCOVERS06

Armrests

Multiple adjustment points are available for the arm brackets. Included is the armrest bracket set with a choice of armrest cushions for left or right hand controller.



Armrest cushion - 32cm/13" LH short **code: XTA80336** RH short **code: XTA80320**



Armrest cushion - 43cm/17" LH Long code: XTA80321 RH Long code: XTA80335



Gel Ovations Arm pad 3.5x14"

code: XAGO314G

Gel Ovations Arm pad 3.5x16"

code: XAGO316G

Gel Ovations Arm pad 2x12"

code: XAGO212G



Armrest bracket set - Dual post flip up (excl armrest cushions) code: XTA80289



Armrest bracket set flip up without stays (excl armrest cushions) code: XTA80302



Stop for armrest brackets (80289 & 80302) code: XTA80743

Real Leather

Gutter Armrest Cushions

The gutter shaped armrests are ideal for users that require extra support of the forearms.



Armrest - Hemi cushion LH long Width: 12cm. Inside width: 7.5cm. Length: 40cm

code: XTA20247



Armrest - Hemi cushion RH short Width: 12cm. Inside width: 7.5cm. Length: 35cm

code: XTA20244



Left hand set-up available on request

Joystick Controllers

To find the best controller for you see the "TA PG-JSM Matrix"



PG JSM LED controller w/o display (no lights) code: XTA95174



PG JSM LED with lights (for chair fitted with LED lights) code: XTA95253



PG JSM Colour screen for hand control (no lights) code: XTA95172



PG JSM Colour screen for hand control with lights (for chair fitted with LED lights) code: XTA95252-1



Can also be fitted for a chair with LED lights and can control TV, DVD, Heat pump etc.

PG JSM Controller large colour display CJSM2-L code: XTA95501

Same but with added function of Bluetooth to connect phone, computer, mouse pointer etc.

PG Controller large colour display CJSM2 -L Bluetooth **code: XTA95524**



All-round joystick Omni (Indent)



Alternative control - Bracket for dropdown s/a (universal) code: XTA80838



Joystick - Swing away bracket RH code: XTA80407 LH code: XTA80408



Joystick – Handle 4" U Flex shift code: XABPPC102B

Joystick knobs




M48 Mushroom code: XTA95154



Ball 38mm code: XTA95155



Cone code: XTA95298



Softball SB70 code: XTA95300



T-Bar code: XTA95301



Stick S80 code: XTA95305



Chin C38 code: XTA95306



Golf ball code: XTA80653



Standard ball code: XTA95153

Alternative Controls



TA hand control code: XTA97237

TA hand control for CJSM2

code: XTA97343

Bracket for TA hand control

code: XTA80527



TA Headcontrol 3 switch head array

code: XTA80681



PB30 Black button 30mm code: XTA91010

PB50 Black button 50mm

code: XTA91033



Gyroset Vigo

Gyroset[™] Vigo is a proportional, wireless headset that enables users to control their powered wheelchairs smoothly with their head movements.

Giving back mobility to patients with: Spinal cord injuries, ALS, Quadriplegia, Muscular dystrophy, Multiple Sclerosis and many more. code: XTA80910



Main Features:

- > The smoothest proportional head control available on the market
- > Wireless Bluetooth technology
- > Connectable to mobile phone
- > Extremely short training period required
- > Customizable settings (all angles and sensitivity) for each individual's needs in the free software
- > Built-in emergency stop functions

Overview:

- > Compatible with PG's OMNI, OMNI 2
- > Super sensitive use it even with the smallest head movements
- > Works with any external switch you prefer, but our Headrest Sensor is included in the kit
- > Re-calibrate neutral position any time
- > Access seating, speed and profile menu with the infrared sensor, adjust each setting easily
- > Activate deactivate drive mode easily




Omni2

The R-net Omni2 is compatible with a whole range of Specialty Input Devices (SIDs) including: joystick types, switch panels, switched and proportional head controls, sip-and-puff and scanners. **code:** XTA95541

- > R-net specialty control Interface
- > Large, high resolution, 3.5" VGA screen
- > Integrated Bluetooth® wireless technology
- > Infra-Red (IR) control
- > Built-in light sensor
- > USB charging socket
- > High-volume horn

- > Dual-channel inputs, analogue and digital
- > Sip-and-Puff transducer
- > On-board configuration programming
- > Electronics sealed to IPx4
- > Backwards compatible 'Connect & Go'

Attendant Controls





Flight Kit LED Lights & USB Charger

Flight Kits

Flight kits are designed to disconnect the batteries easily. This is an requirement for most airlines.

TA iQ Flight kit
IQ code: XTA97285
TA indoor Flight kit

Indoor code: XTA95527











Storage Bags & Hooks Push Handles & Rearview Mirror











Cup & Device Holders



RAM Suction Cup Base with B Series Balll

code: XRAM-B-224-1U




RAM Drink Holder Self-leveling holder keeps your drink upright regardless of the conditions

code: XRAM-B-132-400U XL code: XRAM-B-132-400U-XL




RAM Cellphone Holder code: XRAM-CELL-UN7BU

Device holders available for mobile phone and tablet, multiple sizes available







Tray - large(50cm) L/H code: XTATrayL-LH Tray - large(50cm) R/H code: XTATrayL-RH Tray - medium(45cm) L/H code: XTATrayM-LH
Tray - medium(45cm) R/H code: XTATrayM-RH

Tray - small(40cm) L/H code: XTATrayS-LH
Tray - small(40cm) R/H code: XTATrayS-RH

Customisations of trays are available contact us to discuss your requirements.

Electric Seat Glide Function Seat Pan Systems

Electric Seat Glide Function

Allowing the seat to glide forward at a max range of 4.5" / 11cm.

Giving you the ability to get closer to counters or sit under tables and desks.

Ideal for when the user is driving a vehicle it allows to get closer to the steering wheel.

code: XTA20710







TA iQ Seat spacer plate 1" code: XTASPS1-VAR
TA iQ Seat spacer plate 2" code: XTASPS2-VAR

TA Indoor Seat Spacer plate 1"

code: XTASPS1-VAR

TA Indoor Seat Spacer plate 2"

code: XTASPS2-VAR



Extender for seat plate 38cm (15") (for 40cm chair width, physical plate measurement 32cm) code: XTA80233

Extender for seat plate 45cm (18") (for 45m chair width, physical plate measurement 40.5cm)

code: XTA80216

Extender for seat plate 50cm (20") (for 50cm chair width, physical plate measurement 45.5cm)

code: XTA80215



TA Anterior wedge (Custom) code: XTA80817



Mid plate for all Chairs (needed for 20" and over) code: XTA80502

Seat Belts Back Canes



Pelvic positioning belt. 2 Point padded. Buckle may be opened by one hand.

Medium code: XTA93223 Large code: XTA93224



Pelvic positioning belt. 4 Point padded. Buckle may be opened by one hand.

Medium code: XTA93227 Large code: XTA93228



Back Canes adjustable width - Tall straight 1" cane design

code: XTA80355



Back Canes adjustable width - Contour design (only for Icon backrest)

code: XTA80298

Contour canes are best with a fixed mounted back not needing angle adjustment.

Chair widths - Drive wheel size options

Standard door jambs allow for clearance of 32"/82cm. However some doors are custom made smaller, ie 27.5"/70cm or even narrower. Changing to a narrower wheel can make the difference to a chairs width.

Chair	TAiQ with Std wheels	TAiQ with Narrow wheels	TA Indoor / TA Junior
Chair width	24.75"/63cm	23.5"/60cm	22"/56cm
Tyre height	14" standard	12.5" code: XTA97270	10" standard
Minimum seat pan to floor height	15"/38cm	15"/38cm	15"/38cm
Max client weight	140Kg	140Kg	120Kg

Note: If the doorways are very narrow & tight, consider the TA Indoor or TA Junior that has a width on only 22"/56cm with a very short wheelbase of 27"/69cm, and a tighter turning radius of 18"/45cm.

TA iQ Junior & TA Junior Seating set-up accessories



Headrest system, cushion + bracket - Jnr 5" (12cm) code: XTA80288



Headrest system, cushion + bracket - Med 7" (18cm) code: XTA80337-1



Junior Drop in Seat Guides (Pair) (VAR)

code: XTA94000VAR Custom sizes available on request



Junior Seat adjustable width and depth. (incl extender plate) code: XTA80291J



Extender for seat plate 38cm (15") code: XTA80233



Junior Back canes adjustable width - contour design (only for Icon backrest)

code: XTA80298J



Back Canes adjustable width - Tall straight 1" cane design (Custom modified for Junior)

code: XTA80355



Junior Arm rest cushion RH Short code: XTA80320J LH Long code: XTA80336J



Junior Flip up armrest bracket (set) code: XTA80302J

Junior chairs are compatible with most accessories in this catalogue, including custom colour and print options, please ask!

Vehicle Docking System



The Dahl wheelchair docking system is a safe way to quickly and independently (for some users) secure a power wheelchair in a vehicle.

A plate attached to the underside of the powerchair slides into the Docking System which is mounted to the floor of the vehicle. The shape of the plate ensures easier docking as it guides the chair into a secure and stable position.

With the use of an electric button the dock will be locked. Once locked, the chair will not twist or move while driving.





Docking Plate

The Docking plate fitted under the chair is very flush. It doesn't protrude or catch kerbs or door frames, unlike pin systems.



Interchangeable Seating Positions

The docking station can be used for either wheelchairs or car seats. It is easy to install a connector to the vehicle's original seat and swap the chairs out when taking turns driving the vehicle.



Road Safety

All TA Electric wheelchairs comply with the NZ safety standards (LVVTA 45-60(02) Disability Transportation Systems) and, in particular, international standard ISO 7176-19 (wheeled mobility devices for use as seats in motor vehicles). TA Service have crash tested their chairs to maximise your safety and minimise your injury risk in the event of a collision

Dahl Docking installation kits are available for a wide range of vehicles. Please check the brochure for compatible models on our website or give Morton & Perry a call.

Parts



Plate for Dahl Bracket code: XTA97189



Dahl locking plate code: XD500561



Docking Station code: XD501750


Find the right backrest and cushion for your chair

Visit our website and our:

PRODUCT SELECTION TOOL or COMPARISON CHART

will guide you to the best fit.





Back Systems

Product Selection Tool CLICK HERE

Comparison Chart CLICK HERE

Seating

Product Selection Tool CLICK HERE

Comparison Chart CLICK HERE









For full specifications see our website. Contact our team for a product demonstration. 0800 238 523 | mortonperry.co.nz